Penumbra Embolization Coils: LP™ Coil and swiftPAC™ MRI Safety Information: MR Conditional

Non-clinical testing has demonstrated the coil is MR Conditional. Coil arrays consisting of one or more coil implants can be scanned safely immediately after implantation under the following conditions.

Parameter	Condition
Static Magnetic Field Strength	1.5T or 3T
Maximum Spatial Field Gradient	3000 Gauss/cm (30 T/m) or less
RF Coil Type	Any
Maximum Whole Body SAR	2.0 W/kg
Maximum Head SAR	3.2 W/kg
Scan Duration	2.0 W/kg whole-body average SAR for 60 minutes of continuous RF (a sequence or back to back series/scan without breaks).

Artifact Information

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the Penumbra coil arrays. Therefore, it may be necessary to optimize parameters for the presence of this implant. In non-clinical testing:

- Artifact associated with a Penumbra coil array was assessed using spin echo and gradient echo sequences. The maximum artifact distance beyond the implant was 3 mm for the spin echo sequence and 8 mm for the gradient echo sequence.
- Artifact associated with a Penumbra coil array in the form of an 8 mm diameter sphere was assessed using a clinical MRA sequence. The maximum artifact distance beyond the implant was 2 mm using this sequence.

 Patient's Name

 Treating Physician

 Phone Number

 Hospital information

 Date of treatment or other information

If you have any questions about MRI scans or Penumbra Embolization Coils, contact your physician.

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