INTRODUCTION

This document gives you information about your newly implanted coil (or coils). Talk to your health care provider if you have any questions or concerns that are not answered in this document. This document is not intended to replace a consultation with your healthcare professional.

IMPORTANT INFORMATION ABOUT DEVICE LIFETIME AND POTENTIAL FOLLOW-UP STUDIES

The coil is a permanent implant. This means the device will stay in your body for your entire life. Make sure all of your healthcare team including your primary care physician or specialist are aware of your coil implant.

Your health care provider may order follow-up studies which may include Magnetic Resonance Imaging (MRI). An MRI uses radio waves and a magnetic field to produce pictures of the body's organs and structures. Special considerations for your health care provider and radiology staff are listed under the section "Warning – Magnetic Resonance Imaging (MRI) Safety Information".

APPOINTMENTS WITH YOUR HEALTH CARE PROVIDER

After your coil implantation procedure, your health care provider may schedule follow-up appointments to see you. Consult with your health care provider for a follow-up treatment plan.

WARNING – MAGNETIC RESONANCE IMAGING (MRI) SAFETY INFORMATION

You should talk to your healthcare provider about your coil implant before undergoing any imaging procedure including MRI. Show them your Patient Implant Card. Your coil is MR conditional. This means the coil has been tested for use in MRIs under specific conditions. Special considerations for your health care provider and radiology staff are located at <https://www.penumbrainc.com/de/mri-compatibility/>.

MATERIALS

The Penumbra Embolization Coil implants contain the following materials:

Material Description	CAS	Substance to Implant ratio (%w/w) High: $\geq 30\%$ Medium: $<30\% - \geq 5\%$ Low: $<5\% - \geq 0.1\%$ *Not included: $<0.1\%$.
Platinum-Tungsten Alloy	7440-06-4	High
	7440-33-7	Medium
Polyethylene	9002-88-4	Low
	68187-11-1	Low
	7440-47-3	Low
	91-17-8	Low
Adhesive	2680-03-7	Low
	N/A	Low
	5888-33-5	Low
	73324-00-2	Low
	2399-48-6	Low
	75980-60-8	Low

Material Description	CAS	Substance to Implant ratio (%w/w) High: $\geq 30\%$ Medium: $<30\% - \geq 5\%$ Low: $<5\% - \geq 0.1\%$ *Not included: $<0.1\%$.
Adhesive	24615-84-7	Low
	7328-17-8	Low
	2530-83-8	Low
	79-10-7	Low
	42978-66-5	Low
	818-61-1	Low
Titanium	7440-32-6	Low
Nitinol**	7440-02-0	Low
	7440-32-6	Low
Platinum-Iridium**	7440-06-4	Low
	7439-88-5	Low

*Substances that are not included are considered below the threshold shown for reporting. A maximum implant mass of 30 mg per cm of coil implant length may be used to calculate discrete substance mass.

**Smart Coil Only

The Penumbra Embolization Coils are not made with natural rubber latex. If you are allergic to nickel or any other material listed above, please discuss this with your physician.

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