INTENDED PURPOSE

The INDIGO Aspiration System is designed to remove emboli and thrombi from the vasculature using continuous aspiration.

DEVICE DESCRIPTION

The INDIGO® Aspiration System is comprised of several devices:

- INDIGO CAT RX Aspiration Catheter
- INDIGO Separator™ 4
- INDIGO Aspiration Tubing
- INDIGO Pump Canister/Tubing
- Penumbra Aspiration Pump

The INDIGO CAT RX Aspiration Catheter is a dual lumen rapid exchange catheter that targets aspiration from the Penumbra Aspiration Pump directly to the thrombus, removing it via the INDIGO Aspiration Tubing and depositing it in the INDIGO Pump Canister. The INDIGO Separator 4 may be used, if needed, to clear the lumen of the INDIGO CAT RX Aspiration Catheter should it become blocked with thrombus. The INDIGO CAT RX Aspiration Catheter is introduced through a guide catheter or long introducer sheath into the coronary or peripheral vasculature and guided over a guidewire to the site of the primary occlusion. The INDIGO CAT RX Aspiration Catheter may be provided with a rotating hemostasis valve and a peelable sheath. The INDIGO Separator 4 is provided with an introducer and torque device. The INDIGO CAT RX Aspiration Catheter and INDIGO Separator 4 are visible under fluoroscopy.

The INDIGO CAT RX Aspiration Catheter has a hydrophilic coating of 20 cm on the distal segment of the catheter shaft.

INDICATION FOR USE

The INDIGO CAT RX Aspiration Catheters and INDIGO Separator 4

As part of the INDIGO Aspiration System, the INDIGO CAT RX Aspiration Catheters and INDIGO Separator 4 are indicated for the removal of fresh, soft emboli and thrombi from vessels in the coronary and peripheral vasculature.

The INDIGO Aspiration Tubing

As part of the INDIGO Aspiration System, the INDIGO Sterile Aspiration Tubing is indicated to connect the INDIGO CAT RX Aspiration Catheters to the Penumbra Aspiration Pump.

Penumbra Aspiration Pump

The Penumbra Aspiration Pump is indicated as a vacuum source for Penumbra Aspiration Systems.

INTENDED PATIENT POPULATION

The intended patient population of the INDIGO Aspiration System is patients with fresh, soft emboli or thrombi in coronary and peripheral arteries.

INTENDED USERS

Intended users for this device are physicians who have received appropriate training in interventional techniques.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- The INDIGO Aspiration System should only be used by physicians who have received appropriate training in interventional techniques.
- Do not advance, retract, or use any component of the INDIGO Aspiration System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device or system as a unit. Unrestrained torquing or forced insertion of the catheter or separator against resistance may result in damage to the device or vessel.
- Do not use the INDIGO Aspiration System with a pump other than the Penumbra Aspiration Pump.

PRECAUTIONS

- The device is intended for single use only. Do not resterilize or reuse.
- Resterilization or reuse of the device may compromise the structural integrity of the device and/or lead to device failure which in turn, may result in patient injury, illness or death.
- Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer/distributor.
- Use prior to the "Use By" date.
- Use the INDIGO Aspiration System in conjunction with fluoroscopic visualization.
- Maintain a constant infusion of appropriate flush solution.
- When performing aspiration, ensure that the INDIGO Aspiration Tubing valve is open for only the minimum time needed to remove thrombus. Excessive aspiration or failure to close the INDIGO Aspiration Tubing valve when aspiration is complete is not recommended.
- The INDIGO Separator 4 is not intended for use as a guidewire. If repositioning of the INDIGO CAT RX Aspiration Catheter is necessary during the revascularization procedure, such repositioning should be performed over an appropriate guidewire using standard guidewire techniques.
- Do not use INDIGO Separator 4 to macerate or retrieve thrombus distal to the catheter tip. INDIGO Separator 4 is intended to be used with INDIGO CAT RX Aspiration Catheter to clear the distal end of the catheter lumen should it be blocked with thrombus.
- Do not use automated high-pressure contrast injection equipment with the INDIGO CAT RX Aspiration Catheter because it may damage the device.
- The INDIGO CAT RX Aspiration Catheter, and INDIGO Separator 4 contain a CMR substance (CAS Number: 7440-48-4).

POTENTIAL ADVERSE EVENTS

Possible complications include, but are not limited to, the following:

- allergic reaction and anaphylaxis from contrast media or device material
- acute vessel occlusion
- air embolism
- anemia
- arrhythmia
- arteriovenous fistula
- death
- emboli
- emergent surgery
- foreign body embolization
- hematoma or hemorrhage at access site
- hemorrhage
- infection
- inflammation
- ischemia
- myocardial infarction
- neurological deficits including stroke
- pseudoaneurysm
- renal impairment or acute renal failure from contrast media
- residual thrombus due to inability to completely remove thrombus
- vessel dissection (intimal disruption)
- vessel perforation
- vessel spasm
- vessel thrombosis

*If a serious incident related to the device occurs, contact your Penumbra, Inc. representative and the competent authority in your respective country/region.

PROCEDURE

1. Refer to **Warnings, Precautions, and Potential Adverse Events** prior to use. As each device of the INDIGO Aspiration System is used, remove the device from the packaging, and inspect for damage or for kinks.
2. Prepare INDIGO Aspiration System devices for use by flushing the packaging and device with heparinized saline (includes INDIGO CAT RX Aspiration Catheter and INDIGO Separator 4 only).
3. Prepare a guide catheter or long introducer sheath according to the manufacturer's *Instructions for Use*.
4. Using conventional catheterization techniques under fluoroscopic guidance, place the guide catheter or long introducer sheath into the appropriate vessel that is proximal to the thrombus occlusion site over an appropriate guidewire.

INDIGO CAT RX ASPIRATION CATHETER PREPARATION AND USE

1. Backload an appropriate guidewire through the guidewire lumen and insert the INDIGO CAT RX Aspiration Catheter into a rotating hemostasis valve connected to the hub of a guide catheter. If a guide catheter is not used, insert the INDIGO CAT RX Aspiration Catheter through the valve of the long introducer sheath using a peelable sheath. After inserting the INDIGO CAT RX Aspiration Catheter, remove the peelable sheath from the introducer sheath, and peel from the INDIGO CAT RX Aspiration Catheter shaft.
2. Using fluoroscopy, advance the INDIGO CAT RX Aspiration Catheter into the target vessel over the guidewire. Position the INDIGO CAT RX Aspiration Catheter proximal to the thrombus.

PENUMBRA ASPIRATION PUMP, INDIGO ASPIRATION TUBING, AND INDIGO SEPARATOR 4 PREPARATION AND USE


1. Attach the INDIGO Aspiration Tubing to the Penumbra Aspiration Pump and turn on the Pump (Refer to the Penumbra Aspiration Pump Operation Manual). Ensure the valve on the INDIGO Aspiration Tubing is switched to the OFF position.
2. If needed, select the INDIGO Separator 4 and attach a rotating hemostasis valve to the hub of the INDIGO CAT RX Aspiration Catheter. Using the introducer sheath provided, insert the INDIGO Separator 4 through the rotating hemostasis valve and into the hub of the INDIGO CAT RX Aspiration Catheter. Connect the INDIGO Aspiration Tubing to the side port of the rotating hemostasis valve as shown in **Figure 1**.

Figure 1 – Assembled INDIGO Aspiration System with INDIGO Separator 4

3. Connect the INDIGO Aspiration Tubing to the hub of the catheter.
4. To begin aspiration, turn the INDIGO Aspiration Tubing valve to the ON position, advance the INDIGO CAT RX Aspiration Catheter to gently embed the distal tip into the thrombus, and if needed, advance/retract the INDIGO Separator 4 to assist with aspiration and removal of the thrombus as shown in **Figure 2**.

Figure 2 – Thrombus Aspiration Using the INDIGO Separator 4

5. To stop aspiration, switch the valve on the INDIGO Aspiration Tubing to the OFF position and turn off the Penumbra Aspiration Pump.
6. If used, remove the INDIGO Separator 4.
7. If necessary, additional INDIGO Separator 4 and INDIGO CAT RX Aspiration Catheters may be used to further remove thrombus at the discretion of the physician.
8. Obtain a post-treatment angiogram by injecting contrast media through the guide catheter or long introducer sheath.

CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS

Clinical benefit of the INDIGO Aspiration System is the removal of thrombus from the coronary or peripheral vasculature, to prevent end organ ischemia or infarction. The detailed listing of all known benefits is included in the Summary of Safety and Clinical Performance (SSCP).


SAFETY AND CLINICAL PERFORMANCE

A copy of the Summary of Safety and Clinical Performance can be viewed by searching the device name on the EUDAMED website: <https://ec.europa.eu/tools/eudamed/eudamed>
Basic UDI DI: 081454801CRXES

STORAGE

- Store in a cool, dry place.
- Dispose of device in accordance with hospital, administrative, and/or local government policy.

SYMBOLS GLOSSARY

	Attention, see instructions for use
	Prescription only – US Federal Law restricts this device to use by or on the order of a physician
	Nonpyrogenic
	Sterilized using ethylene oxide. Single sterile barrier system with protective packaging outside
	Medical Device
 Manufacturer	 EU Representative
 Lot number	 Do not reuse
 Use by	 Catalogue number
 Do not resterilize	 Not made with natural rubber latex
 Maximum Guidewire Diameter	 6F Guide Catheter Compatible
 Do not use if package is damaged.	 Contains Hazardous Substances CAS Number: 7440-48-4

WARRANTY

Penumbra, Inc. (Penumbra) warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning, and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Penumbra's control directly affect the device and the results obtained from its use. Penumbra's obligation under this warranty is limited to the repair or replacement of this device and Penumbra shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. Penumbra neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. Penumbra assumes no liability with respect to devices reused, reprocessed, or reesterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such device.

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