

## Swift Coil System MRI Information: MR Conditional

### MRI Safety Information

Non-clinical testing has demonstrated the Penumbra Swift Coil Implant is MR Conditional. Coil arrays consisting of one or more Penumbra Swift Coil Implants can be scanned safely immediately after implantation under the following conditions.

Parameter	Condition
Static Magnetic Field Strength	1.5T or 3T
Maximum Spatial Field Gradient	3000 Gauss/cm (30 T/m) or less
RF Coil Type	Any
Maximum Whole Body SAR	2.0 W/kg
Maximum Head SAR	3.2 W/kg
Scan Duration	2.0 W/kg whole-body average SAR for 60 minutes of continuous RF (a sequence or back to back series/scan without breaks).

**Swift Coil System MRI Information: See back of card for more information.**

---

Patient's Name

---

Treating Physician

Phone Number

---

Hospital Information

---

Date of treatment or other information

If you have any questions about MRI scans or Swift Coil System, contact your physician.

Penumbra



[penumbrainc.com](http://penumbrainc.com)