**DEVICE DESCRIPTION**

Penumbra swiftPAC™ Coil consists of a platinum embolization coil attached to a composite detachment pusher and Detachment Handle. The Detachment Handle is packaged separately.

**INTENDED PURPOSE**

Penumbra swiftPAC Coil is intended to endovascularly obstruct or occlude blood flow in aneurysms or other vascular abnormalities.

**INDICATION FOR USE**

Penumbra swiftPAC Coil is indicated for the embolization of:

- Intracranial aneurysms
- Other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae
- Arterial and venous embolizations in the peripheral vasculature

**CONTRAINDICATIONS**

There are no known contraindications.

**INTENDED PATIENT POPULATION**

The intended patient population includes patients harboring aneurysms and other vascular malformations.

**INTENDED USERS**

Intended users for this device are physicians who have received appropriate training in interventional techniques.

**DEVICE COMPATIBILITY**

The device requires a delivery microcatheter with a minimum inner diameter of 0.0165 in. [0.42 mm].

**RECOMMENDED ACCESSORY PRODUCTS**

1. Femoral sheath
2. Guide catheter compatible with microcatheter and any adjunct devices
3. Microcatheter with a compatible inner diameter
4. Guide wire compatible with microcatheter
5. Rotating Hemostasis Valves (RHVs)
6. Accessories for continuous flush set-up
7. Extra Detachment Handle (recommended)

**WARNINGS**

- This device should only be used by physicians who have received appropriate training in interventional techniques.
- Do not use kinked or damaged devices. Do not use opened or damaged packages. Use of damaged coils may affect coil delivery and/or stability inside the aneurysm or vessel. Return damaged devices and packaging to the manufacturer/distributor.
- Do not advance or withdraw the device against resistance without careful assessment of the cause using fluoroscopy. Moving or torquing the device against resistance may result in damage to the vessel or device.
- Remove the coil if unusual friction is encountered. If resistance is encountered when withdrawing the coil, withdraw the microcatheter until the resistance subsides. If unusual friction is encountered on subsequent coils, remove the microcatheter and carefully examine for damage.
- Do not rotate the delivery pusher during use. Rotating the delivery pusher may result in premature detachment, which could lead to coil damage, incorrect positioning, or vessel damage.
- Verify repeatedly that the microcatheter is not under stress before coil detachment. Stored forces in the microcatheter could cause the tip to move during detachment, which could lead to lesion rupture.
- Advancing the delivery pusher beyond the microcatheter tip could lead to lesion rupture.

**PRECAUTIONS**

- The device is intended for single use only. Do not resterilize or reuse. Resterilization and/or reuse may compromise the structural integrity of the device or increase the risk of contamination or infection leading to device failure and/or cross-infection and potential patient injury, illness or death.
- Use prior to the "Use By" date.
- Use device in conjunction with fluoroscopic guidance.
- Do not advance or retract the device against resistance without careful assessment of the cause using fluoroscopy. If device cannot be advanced or retracted, withdraw the device as a unit with the microcatheter.
- As in all fluoroscopy procedures, consider all necessary precautions to limit patient radiation exposure by using sufficient shielding, reducing fluoroscopy times and modifying radiation technical factors whenever possible.
- Maintain a constant infusion of an appropriate flush solution. Lack of an RHV and/or continuous infusion of appropriate flush solution into the microcatheter may result in difficulty transferring the coil from the introducer sheath and/or difficulty advancing the coil through the microcatheter, which may lead to kinking or coil detachment.
- The device may create local field inhomogeneity and susceptibility artifacts during magnetic resonance angiography (MRA), which may degrade the diagnostic quality to assess effective lesion treatment.
- Penumbra Embolization Coil System detachment pusher contains a CMR substance (CAS Number: 7440-48-4).

**POTENTIAL ADVERSE EVENTS**

Potential complications include but are not limited to:

- |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>• acute vessel occlusion</li> <li>• air embolism</li> <li>• allergic reaction and anaphylaxis from contrast media or device materials</li> <li>• aneurysm rupture</li> <li>• arteriovenous fistula</li> <li>• coil herniation into parent vessel</li> <li>• death</li> <li>• emboli (includes foreign body embolization)</li> <li>• pseudoaneurysm</li> <li>• hematoma or hemorrhage at access site</li> <li>• incomplete aneurysm occlusion requiring retreatment</li> <li>• infection</li> <li>• inflammation</li> <li>• intracranial hemorrhage</li> </ul> | <ul style="list-style-type: none"> <li>• ischemia</li> <li>• myocardial infarction</li> <li>• necrosis</li> <li>• neurological deficits including stroke</li> <li>• thromboembolic events</li> <li>• post-embolization syndrome</li> <li>• recanalization</li> <li>• renal impairment or acute renal failure from contrast media</li> <li>• vessel spasm, thrombosis, dissection (intimal dissection) or perforation</li> <li>• radiation exposure that may lead to cataracts, skin reddening, burns, alopecia, or neoplasia from x-ray exposure</li> </ul> |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

If a serious incident related to the device occurs, contact your Penumbra, Inc. representative and the competent authority in your respective country/region.

**DEVICE PREPARATION AND USE**

Refer to Warnings, Precautions and Potential Adverse Events prior to use.

**CONTINUOUS FLUSH SET UP**

To maintain optimal performance and reduce the risk of thromboembolic events, it is critical that a continuous infusion of appropriate flush solution is maintained into the femoral sheath, guide catheter, and microcatheter.

- Attach a RHV to the guide catheter hub.
- Connect a three way stopcock to the RHV sidearm. Attach a continuous flush line to the stopcock.
- Attach a RHV to the microcatheter hub.
- Connect a one way stop cock to the RHV sidearm. Attach a continuous flush line to the stopcock.
- Insert the microcatheter into the RHV connected to the guide catheter hub.
- Adjust continuous flush pressure to maintain appropriate infusion of solution during procedure.
- Check all fittings to ensure no air is introduced into system.

**CATHETERIZATION**

- Using conventional catheterization techniques under fluoroscopic guidance, access the target vessel with the guide catheter. The guide catheter should have an inner diameter large enough to allow for sufficient contrast injection during the procedure.
- Carefully catheterize the lesion with the microcatheter according to the manufacturer's instructions. After the microcatheter is in position, remove the guide wire.

**COIL SELECTION**

- Perform fluoroscopic road mapping.
- Measure the size of the lesion to be treated.
- Select an appropriate size and configuration for the initial coil.
- The diameter of the first coil should generally be larger than the width of the aneurysm neck or the propensity for coils to migrate may increase.
- In order to choose the optimum coil for any given lesion, examine pretreatment angiograms. The appropriate coil size should be chosen based upon angiographic assessment of the dimensions, overall morphology, and rupture status of the lesion and surrounding vasculature.
- Subsequent coils may be selected, depending on preceding coil deployment and the remaining lesion volume to be filled. Coils will typically be of decreasing size until the physician determines that the lesion has been successfully treated.


**DIRECTIONS FOR USE**

1. Slowly remove the coil and introducer sheath from dispenser hoop by gripping both the pusher and introducer sheath. Use caution to avoid damaging the delivery pusher. Carefully examine the coil delivery pusher for damage or irregularities.
2. Carefully advance the coil out of the introducer sheath into heparinized saline. Examine the coil for kinks or other irregularities. If an irregularity is found, replace with a new coil.
3. With the distal (tapered) end of the introducer sheath pointed downward and submerged in saline, carefully withdraw the coil fully into the introducer sheath.
4. Introduce the distal end of the introducer sheath through the RHV into the hub of the microcatheter until the introducer sheath is firmly seated in the hub taper. Tighten the RHV around the introducer sheath to prevent retrograde flow of blood. Do not over-tighten. Excessive tightening could damage the coil or delivery pusher.
5. Transfer the coil into the microcatheter by gently advancing the delivery pusher in a smooth continuous manner. Once the flexible portion of the delivery pusher has entered the microcatheter, loosen the RHV and remove the introducer sheath over the delivery pusher's proximal end.

**NOTE: To resheath the coil, backload the delivery pusher into the distal (tapered) end of the introducer sheath. With the distal end of the introducer sheath pointed downward, carefully withdraw the coil fully into the introducer sheath.**


6. Tighten the RHV on the delivery pusher to prevent retrograde blood flow while still allowing advancement of the delivery pusher.
7. Visually verify that the flushing solution is infusing normally.
8. When using a 150 cm or longer length microcatheter, advance the delivery pusher until the coil fluoroscopy marker reaches the RHV. Fluoroscopic guidance should be initiated at this time. See **Figure 1**. If using microcatheter of a different length this entire step should be done under fluoroscopic guidance.

**Figure 1: Coil Fluoroscopy Marker Reaches the RHV**

9. Advance the coil under fluoroscopic guidance to the desired site. If coil placement is not satisfactory, slowly withdraw by pulling the delivery pusher proximally. The coil can then be deployed again. If the coil size is inappropriate, remove and replace with an appropriately sized coil.
10. Continue advancing the coil until the radiopaque marker on the delivery pusher aligns just beyond the microcatheter proximal radiopaque marker to form a "T". At this point, the proximal end of the coil implant just exits the microcatheter tip. Tighten the RHV on the delivery pusher to prevent further movement. See **Figure 2**.

**Figure 2: Microcatheter Proximal Marker Forms A "T" With Delivery Pusher Marker**

11. Remove the Detachment Handle from its packaging and place it within the sterile field.
- NOTE: Verify that the Detachment Handle is labeled for use with corresponding family of coils.**
12. Reconfirm the delivery pusher radiopaque marker forms a “T” with the microcatheter proximal marker.
- NOTE: Coil detachment is not dependent on exact alignment. However, marker alignment will ensure that the proximal end of the coil implant just exits the microcatheter tip.**
13. Hold the delivery pusher distal to the black alignment zone. Advance the Detachment Handle over the proximal end of the delivery pusher until the black alignment zone is fully concealed within the Detachment Handle funnel and the delivery pusher comes to a stop inside the Detachment Handle. See **Figure 3**.

Figure 3: Delivery Pusher Aligns With Detachment Handle



Figure 4: Detachment Handle

14. To detach the coil, use thumb (**Figure 4**) to retract the Detachment Handle slider until it clicks.
- NOTE: After the slider clicks, release the slider to allow it to return to its original position.**
- NOTE: The detachment mechanism will only activate upon hearing the click, not during the earlier phase of slider retraction.**
15. Gently pull back on the Detachment Handle to remove it from the delivery pusher. Gently pull back the delivery pusher under fluoroscopic guidance to ensure the coil has detached.
16. If the coil is still attached, repeat steps 13 through 15, or carefully withdraw the coil and deploy another coil.
17. Repeat for all additional coils.
18. Once the procedure is complete, discard the Detachment Handle.

MRI SAFETY INFORMATION

|                                                                                                                                                                                                            |                                                                                                                          |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------|
| Non-clinical testing has demonstrated the coil is MR Conditional. Coil arrays consisting of one or more coil implants can be scanned safely immediately after implantation under the following conditions. |                                                                                                                          |
| Parameter                                                                                                                                                                                                  | Condition                                                                                                                |
| Static Magnetic Field Strength                                                                                                                                                                             | 1.5T or 3T                                                                                                               |
| Maximum Spatial Field Gradient                                                                                                                                                                             | 3000 Gauss/cm (30 T/m) or less                                                                                           |
| RF Coil Type                                                                                                                                                                                               | Any                                                                                                                      |
| Maximum Whole Body SAR                                                                                                                                                                                     | 2.0 W/kg                                                                                                                 |
| Maximum Head SAR                                                                                                                                                                                           | 3.2 W/kg                                                                                                                 |
| Scan Duration                                                                                                                                                                                              | 2.0 W/kg whole-body average SAR for 60 minutes of continuous RF (a sequence or back to back series/scan without breaks). |

ARTIFACT INFORMATION

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the Penumbra coil arrays. Therefore, it may be necessary to optimize parameters for the presence of this implant. In non-clinical testing:

- Artifact associated with a Penumbra coil array was assessed using spin echo and gradient echo sequences. The maximum artifact distance beyond the implant was 3 mm for the spin echo sequence and 8 mm for the gradient echo sequence.
- Artifact associated with a Penumbra coil array in the form of an 8 mm diameter sphere was assessed using a clinical MRA sequence. The maximum artifact distance beyond the implant was 2 mm using this sequence.

CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS

Penumbra Embolization Coil System endovascularly obstructs or occludes blood flow in neurovascular or peripheral blood vessels to prevent aneurysm rupture or bleeding from vascular abnormalities. The detailed listing of all known benefits is included in the Summary of Safety and Clinical Performance (SSCP).

INTENDED DEVICE LIFETIME

The intended lifetime of the coil implant is the lifetime of the patient.

PACKAGING AND STORAGE

The device is packaged in a protective plastic dispenser hoop, pouch, and box. The device will remain sterile unless the pouch is opened, damaged, or the “Use By” date has passed. The Detachment Handle is packaged separately in a protective pouch and box. The Detachment Handle will remain sterile unless the pouch is opened, damaged, or the “Use By” date has passed.

MATERIALS

The following materials are present in the swiftPAC Coil implants (see **Table 1**) and have a substance to material ratio of >0.1% (%w/w). Where applicable, CAS numbers are provided.

| Table 1: Coil Implant Materials to Which the Patient may be Exposed |            |                                                                                                                   |
|---------------------------------------------------------------------|------------|-------------------------------------------------------------------------------------------------------------------|
| Material Description                                                | CAS        | Substance to Implant ratio (%w/w)<br>High: ≥30%<br>Medium: <30% - ≥5%<br>Low: <5% - ≥0.1%<br>*Not included: <0.1% |
| Platinum-Tungsten Alloy                                             | 7440-06-4  | High                                                                                                              |
|                                                                     | 7440-33-7  | Medium                                                                                                            |
| Polyethylene                                                        | 9002-88-4  | Low                                                                                                               |
|                                                                     | 68187-11-1 | Low                                                                                                               |
|                                                                     | 7440-47-3  | Low                                                                                                               |
|                                                                     | 91-17-8    | Low                                                                                                               |
| Adhesive                                                            | 2680-03-7  | Low                                                                                                               |
|                                                                     | N/A        | Low                                                                                                               |
|                                                                     | 5888-33-5  | Low                                                                                                               |
|                                                                     | 73324-00-2 | Low                                                                                                               |
|                                                                     | 2399-48-6  | Low                                                                                                               |
|                                                                     | 75980-60-8 | Low                                                                                                               |
|                                                                     | 24615-84-7 | Low                                                                                                               |
|                                                                     | 7328-17-8  | Low                                                                                                               |
|                                                                     | 2530-83-8  | Low                                                                                                               |
|                                                                     | 79-10-7    | Low                                                                                                               |
| Titanium                                                            | 42978-66-5 | Low                                                                                                               |
|                                                                     | 818-61-1   | Low                                                                                                               |
|                                                                     | 7440-32-6  | Low                                                                                                               |

\*Substances that are not included are considered below the threshold shown for reporting unless disclosed per MDR Section 10.4.5 (≥0.1% w/w) or other applicable regulations or standards. A maximum implant mass of 30 mg per cm of coil implant length may be used to calculate discrete substance mass.

The Coils, Detachment Handles, and packaging are not made with natural rubber latex.

SAFETY AND CLINICAL PERFORMANCE

A copy of the Summary of Safety and Clinical Performance can be viewed by searching the device name on the EUDAMED website: <https://webgate.ec.europa.eu/eudamed>.

BASIC UDI-DI:

A complete list of the basic UDI-DI assigned to the components of the Penumbra Embolization Coil System is provided in **Table 2**.

| Table 2: Penumbra Embolization – List of Basic UDI-DI |                                                    |                                                       |
|-------------------------------------------------------|----------------------------------------------------|-------------------------------------------------------|
| System                                                | Basic UDI-DI Groups                                | Basic UDI-DI (GS1 Prefix + Description + Check Digit) |
| Penumbra Embolization Coil System                     | Penumbra Embolization Coils and Detachment Handles | 081454801PCOEU                                        |

STORAGE AND DISPOSAL

- Store in a cool dry place.
- Dispose of device in accordance with hospital, administrative, and/or local government policy.

SYMBOLS GLOSSARY

|                                                                                   |                                                                                                   |                                                                                                   |                                                        |
|-----------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------|--------------------------------------------------------|
|  |                                                                                                   | Prescription only – US Federal Law restricts this device to use by or on the order of a physician |                                                        |
|  |                                                                                                   | Date of Manufacture                                                                               |                                                        |
|  |                                                                                                   | 10 Microcatheter Compatible                                                                       |                                                        |
|   | Do not use if package is damaged                                                                  |                  | Contains Hazardous Substances<br>CAS Number: 7440-48-4 |
|   | Attention, see Instructions for Use                                                               |                  | Nonpyrogenic                                           |
|   | Sterilized using ethylene oxide. Single sterile barrier system with protective packaging outside. |                  | Medical Device                                         |
|   | Manufacturer                                                                                      |                  | EU Representative                                      |
|   | Lot number                                                                                        |                  | Do not reuse                                           |
|   | Use by                                                                                            |                  | Catalogue number                                       |
|   | Do not resterilize                                                                                |                  | Not made with natural rubber latex                     |

WARRANTY

Penumbra, Inc. (Penumbra) warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning, and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Penumbra's control directly affect the device and the results obtained from its use. Penumbra's obligation under this warranty is limited to the repair or replacement of this device and Penumbra shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. Penumbra neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. Penumbra assumes no liability with respect to devices reused, reprocessed, or re-sterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such devices.

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