

High Rates of Technical Success Indigo System for Acute Limb Ischemia

STRIDE Study


Objective: To collect safety and performance data on the Indigo Aspiration System in a patient population

with lower extremity acute limb ischemia (LE-ALI)

Prospective | Single-arm | Multicenter, 16 Sites | 119 Patients Enrolled Design:











| Baseline Characteristics | % (n/N) or Mean ± SD Median [IQR] | |
|----------------------------|--|--|
| Target Thrombus Length, mm | 125.7 ± 124.7 (N=113) | |
| Rutherford IIa and IIb | 89.1 % (106/119) | |
| Tandem Lesion | 18.5 % (22/119) | |

Patients with firstline use of the **Indigo** Aspiration System had excellent limb salvage rates and low peri-procedural complications.

In this study, Indigo Aspiration System provided a safe and effective clot removal option for LE-ALI patients.

- a. STRIDE study was not a randomized or head-to-head study. Please refer to specific publications to review source for detailed patient and data collection methods for open surgical revascularization. b. Composite limb salvage rate at 30 days calculated and data on file at Penumbra, Inc.

- 1. Maldonado TS, Powell A, Wendorff H, et al. Safety and efficacy of mechanical aspiration thrombectomy for patients with acute lower extremity ischemia. J Vasc Surg. 2024;79(3):584–592.e5. doi:10.1016/j.jvs.2023.10.062. 2. Veenstra EB, van der Laan MJ, Zeebregts CJ, et al. A systematic review and meta-analysis of endovascular and surgical revascularization techniques in acute limb ischemia. J Vasc Surg. 2020 Feb;71(2):654–668.e3. doi:10.1016/j.jvs.2019.05.031.
- 3. Taha AG, Byrne RM, Avgerinos ED, et al. Comparative effectiveness of endovascular versus surgical revascularization for acute lower extremity ischemia. J Vasc Surg. 2015 Jan;61(1):147–154. doi:10.1016/j.jvs.2014.06.109. The clinical results presented herein are for informational purposes only, and may not be predictive for all patients. Individual results may vary depending on patient-specific attributes and other factors. Please refer to the Instructions for Use (IFU) for complete product indications, contraindications, warnings, precautions, potential adverse events, and detailed instructions for use



Historical Surgical Outcomes vs. STRIDE Study Data

| Outcome | Open Surgery | STRIDE ⁶ |
|---------------------------------|---------------------------|-------------------------------------|
| Target Limb Salvage at 30 days | 83.1% ^{2,b} | 98.2% ¹ (109/111) |
| Target Limb Salvage at 365 days | 77.3% ^{2,b} | 88.5 % (77/87) |
| Patency at 30 days | 78.6%' | 89.4 %¹ (101/113) |
| Mortality at 30 days | 13.2%³ | 3.4 %¹ (4/119) |
| Mortality at 365 days | 33.8%³ | 12.0 % ⁵ (12/100) |
| Major Bleeding ^e | 21.0%° | 4.2% ¹ (5/119) |
| Access Site Infection Rate | 14.0% ⁷ | 0.8% ^d |

a. STRIDE study was not a randomized or head-to-head study. Please refer to specific publications to review source for detailed patient and data collection methods for open surgical revascularization.

Penumbra, Inc. USA One Penumbra Place Alameda, CA 94502 USA 1.888.272.4606 T 1.510.748.3200 F 1.510.748.3232

info@penumbrainc.com

Penumbra Europe GmbH Am Borsigturm 44 13507 Berlin Germany T +49 30 2005 676-0 F +49 30 2005 676-10 de-order@penumbrainc.com de-info@penumbrainc.com Penumbra Neuro Australia Pty Ltd 55 Kirby Street Rydalmere NSW 2116 Australia T +61-1300 817 025 F +61-1300 817 026 order.anz@penumbrainc.com

Penumbra Latin America Distribuidora de Equipamentos e Produtos Médicos Ltda Av. Brigadeiro Faria Lima 1336 Cj 82, CEP 01451-001 T +55 11 2883-5825 order.la@penumbrainc.com



For the complete Penumbra IFU Summary Statements, please scan QR code or visit: peninc.info/risk

The clinical results presented herein are for informational purposes only, and may not be predictive for all patients. Individual results may vary depending on patient-specific attributes and other factors.

Product availability varies by country. Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician. Prior to use, please refer to the instructions for Use (IFU) for complete product indications, contraindications, warnings, precautions, potential adverse events, and detailed instructions for use. Please contact your local Penumbra representative for more information.

Copyright ©2023–2024 Penumbra, Inc. All rights reserved. The Penumbra P logos, Indigo, CAVT, Lightning, and Penumbra ENGINE are registered trademarks or trademarks of Penumbra, Inc. in the USA and other countries. 25163 Rev D 09/24 UNI



<sup>b. Composite limb salvage rate at 30 days calculated and data on file at Penumbra, Inc.
c. Major bleeding definitions may vary across studies. Please refer to specific publications for details.
d. Data on file at Penumbra, Inc. based on 30-day outcomes. The clinical results presented herein are for informational purposes.</sup>

^{1.} Maldonado TS, Powell A, Wendorff H, et al. Safety and efficacy of mechanical aspiration thrombectomy for patients with acute lower extremity ischemia. J Vasc Surg. 2024;79(3):584–592.e5. doi:10.1016/j.jvs.2023.10.062. 2. Veenstra EB, van der Laan MJ, Zeebregts CJ, et al. A systematic review and meta-analysis of endovascular and surgical revascularization techniques in acute limb ischemia. J Vasc Surg. 2020 Feb;71(2):654–668.e3. doi:10.1016/j.jvs.2019.05.031.

^{3.} Taha AG, Byrne RM, Avgerinos ED, et al. Comparative effectiveness of endovascular versus surgical revascularization for acute lower extremity ischemia. J Vasc Surg. 2015 Jan;61(1):147–154. doi:10.1016/j.jvs.2014.06.109. 4. Grip O, Wanhainen A, Michaëlsson K, Lindhagen L, Björck M. Open or endovascular revascularization in the treatment of acute lower limb ischaemia. Br J Surg. 2018 Nov;105(12):1598–1606. doi:10.1002/bjs.10954. 5. Maldonado TS, Powell A, Wendorff H, for the STRIDE Study group. One-year limb salvage and quality of life following mechanical aspiration thrombectomy in patients with acute lower extremity ischemia. J Vasc Surg. 2024.

uoi. 10.1016/j.jvs.2024-0.00-03.

K. Kolte D, Kennedy KF, Shishehbor MH, et al. Endovascular versus surgical revascularization for acute limb ischemia: a propensity-score matched analysis. Circ Cardiovasc Interv. 2020;13(1):e008150.

R. Bakshi SC, Fobare A, Benarroch-Gampel J, Teodorescu V, Rajani RR. Lower socioeconomic status is associated with groin wound complications after revascularization for peripheral artery disease. Ann Vasc Surg. 2020 Jan;62:76-82. doi:10.1016/j.avsg.2019.05.010. Epub 2019 Jun 12.