

CAVT™ · The Most Advanced Treatment for PE



STORM-PE RCT

**CAVT · Proven by
Level 1 Evidence**

In Collaboration with
The PERT Consortium™




CAVT · The Most Advanced Treatment for PE

CAVT · Proven by Level 1 Evidence



CAVT shows superior efficacy over anticoagulation alone^a




Objective: Evaluate the efficacy and assess the safety of treating acute, intermediate-high risk pulmonary embolism with anticoagulation plus CAVT Computer Assisted Vacuum Thrombectomy with the Indigo™ Aspiration System Lightning Flash™ versus anticoagulation alone

Design: 100 patients randomized 1:1 to CAVT plus anticoagulation (CAVT arm) or anticoagulation alone (AC arm); 22 U.S. and international sites


Primary Endpoints (Δ RV/LV)

The CAVT arm demonstrated a significantly greater reduction in right heart strain compared to the AC arm

RV/LV Reduction Baseline to 48 Hours



48 Hour RV/LV Reduction > 0.2




a. Efficacy was predefined as the difference between treatment arms in change of RV/LV ratio from baseline to 48 hrs. STORM-PE demonstrated superiority to anticoagulation utilizing Lightning Flash 1.0 and 2.0.

Safety

The CAVT arm safety rate was comparable to the AC arm, with numerically fewer Major Adverse Events (MAE)^b in the CAVT arm^c

Composite MAE ≤ 7 Days (% , n)



Procedural Data · CAVT Arm

Time ^e (median)	mPAP Reduction (mean)	Technical Success ^f	Device Related Transfusion
Device: 25 min Procedure: 56 min	27.3% 8.2 mmHg	100% (47/47)	0% (0/47)

b. MAEs within 7 days included: clinical deterioration necessitating rescue therapy, PE-related mortality, symptomatic recurrent PE, and major bleeding. c. STORM-PE was not powered to detect differences in safety. d. Death, clinical deterioration, and major bleeding all occurred in the same CAVT patient. Type 3a was not considered major bleeding if it was related to an expected decrease in hemoglobin level due to fluid administration and if transfusion was less than 2 units. e. IQR for Device and Procedure time were [15.0, 41.0] and [42.0, 96.0], respectively. f. Technical Success was defined as ability of catheter to access clot and perform aspiration.

DATA NOT YET PUBLISHED. Presented by Lookstein, R. STORM-PE a prospective, multicenter, randomized controlled trial evaluating anticoagulation alone vs. anticoagulation plus mechanical aspiration with the Indigo aspiration system for the treatment of intermediate high risk acute pulmonary embolism. Presented at: TCT (Transcatheter Cardiovascular Therapeutics); October 26, 2025; San Francisco, CA, USA.

Baseline Clinical Parameters^g

Both groups were well matched across key baseline measures

	CAVT Arm (N=47)	AC Arm (N=53)
NEWS2 ^h	3.5 ± 1.95	4.1 ± 2.07
Heart Rate (bpm)	93.2 ± 17.36	98.2 ± 15.87
Oxygen Saturation (%)	96.0 ± 2.59	95.4 ± 2.44
RV/LV Ratio	1.63 ± 0.36	1.56 ± 0.35
RMMS ^h	27.3 ± 3.89	26.1 ± 5.51

Baseline Functional Assessments^g

Borg Dyspnea Scale	4.6 ± 2.90	4.5 ± 2.85
mMRC ^h ≥ 1	43 (93.5%)	51 (96.2%)
PVFS ^h	3.0 ± 0.97	2.9 ± 1.01

g. Paired data at 48 h, paired data for CAVT ranged from N=45–46 and for AC=52–53. Data presented as mean ± SD.
h. NEWS2 = National Early Warning Score 2; RMMS = Refined Modified Miller Score; mMRC = Modified Medical Research Council Dyspnea Scale; PVFS = Post-VTE (Venous Thromboembolism) Functional Status Scale.

More to Come at
VIVA 2025
Las Vegas, NV, USA



Learn More About
STORM-PE RCT

For the complete Penumbra
IFU Summary Statements
and more, scan QR code or
visit: peninc.info/storm-rect



Penumbra, Inc. USA
One Penumbra Place
Alameda, CA 94502
USA
1.888.272.4606
T 1.510.748.3200
F 1.510.748.3232
order@penumbrainc.com
info@penumbrainc.com
penumbrainc.com

Penumbra Europe GmbH
Am Borsigturm 44
13507 Berlin
Germany
T +49 30 2005 676-0
F +49 30 2005 676-10
de-order@penumbrainc.com
de-info@penumbrainc.com

Penumbra Neuro Australia Pty Ltd
55 Kirby Street
Rydalmere NSW 2116
Australia
T +61-1300 817 025
F +61-1300 817 026
order.anz@penumbrainc.com

Penumbra Latin America
Distribuidora de Equipamentos
e Produtos Médicos Ltda
Av. Brigadeiro Faria Lima 1336
Cj 82, CEP 01451-001
São Paulo/SP, Brazil
T +55 11 2883-5825
order.la@penumbrainc.com

The clinical results presented herein are for informational purposes only, and may not be predictive for all patients. Individual results may vary depending on patient-specific attributes and other factors. Product availability varies by country. Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician. Prior to use, please refer to the Instructions for Use (IFU) for complete product indications, contraindications, warnings, precautions, potential adverse events, and detailed instructions for use. Please contact your local Penumbra representative for more information.

Copyright ©2025 Penumbra, Inc. All rights reserved. The Penumbra P logos, the STORM-PE logo, CAVT, Indigo, and Lightning Flash are registered trademarks or trademarks of Penumbra, Inc. in the USA and other countries. All other trademarks are the property of their respective owners. 33137, Rev. A 10/25 UNI