


Select™ Catheter Tip Shapes



Available Lengths

5 F	–	130 cm	120 cm 130 cm	120 cm 130 cm
6 F	125 cm	125 cm	105 cm 125 cm	105 cm 125 cm



6 F Select through
Neuron™ MAX



5 F Select through
BENCHMARK™


Seamless transition zone for atraumatic advancement

Select™ Catheter Construction

Penumbra Select Catheters


- Stainless braided shaft with radiopaque polymer steam-shaped tip
- Tapered profile of each Select catheter optimized for seamless transition and torque response
- Penumbra Select catheters do not have hydrophilic coating

5 F Select | 120/130 cm



- Designed to deliver soft-tipped BENCHMARK™ 071 and Neuron™ 070 Intracranial Access Catheters
- Not designed to deliver 6 F Long Sheath

6 F Select | 105/125 cm



- Designed to deliver Neuron MAX Long Sheath
- Maximum OD .087"
- Not designed for use with smaller ID catheters

Neuron Intracranial Access System — Indication for Use
The Neuron Intracranial Access System is indicated for the introduction of interventional devices into the peripheral, coronary, and neuro vasculature.

Contraindications
There are no known contraindications.

Warnings
The Neuron Intracranial Access System should only be used by physicians who have received appropriate training in interventional techniques.

- Precautions**
- The device is intended for single use only. Do not resterilize or reuse.
 - Resterilization and/or reuse may result in ineffective catheter coating lubrication, which may result in high friction and the inability to access the target vasculature location; and/or may compromise the structural integrity of the device.
 - Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer/distributor.
 - Use prior to the "Use By" date.
 - Use the Neuron Intracranial Access System in conjunction with fluoroscopic visualization.
 - Do not advance or withdraw the Neuron Intracranial Access System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device. Unrestrained moving or torquing the device against resistance may result in damage to the vessel or device.
 - Maintain a constant infusion of an appropriate flush solution.
 - If flow through the device becomes restricted, do not attempt to clear the lumen by infusion. Remove and replace the device.

Potential Adverse Events
Possible complications include, but are not limited to, the following: acute occlusion; air embolism; death; distal embolization; emboli; false aneurysm formation; hematoma or hemorrhage at puncture site; infection; intracranial hemorrhage; ischemia; neurological deficits including stroke; vessel spasm, thrombosis, dissection, or perforation.

Neuron MAX System — Indication for Use
The Neuron MAX System is indicated for the introduction of interventional devices into the peripheral, coronary, and neuro vasculature.

Contraindications
There are no known contraindications.

Warnings
The Neuron MAX System should only be used by physicians who have received appropriate training in interventional techniques.

- Precautions**
- The device is intended for single use only. Do not resterilize or reuse.
 - Resterilization and/or reuse may result in ineffective catheter coating lubrication, which may result in high friction and the inability to access the target vasculature location; and/or may compromise the structural integrity of the device.
 - Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer/distributor.
 - Use prior to the "Use By" date.
 - Use the Neuron MAX System in conjunction with fluoroscopic visualization.
 - Do not advance or withdraw the Neuron MAX System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device. Unrestrained moving or torquing the device against resistance may result in damage to the vessel or device.
 - Maintain a constant infusion of an appropriate flush solution.
 - If flow through the device becomes restricted, do not attempt to clear the lumen by infusion. Remove and replace the device.

Potential Adverse Events
Possible complications include, but are not limited to, the following: acute occlusion; air embolism; death; distal embolization; emboli; false aneurysm formation; hematoma or hemorrhage at puncture site; infection; intracranial hemorrhage; ischemia; neurological deficits including stroke; vessel spasm, thrombosis, dissection, or perforation.

BENCHMARK Intracranial Access System — Indication for Use
The BENCHMARK Intracranial Access System is indicated for the introduction of interventional devices into the peripheral, coronary, and neuro vasculature.

Contraindications
There are no known contraindications.

Warnings
The BENCHMARK Intracranial Access System should only be used by physicians who have received appropriate training in interventional techniques.

- Precautions**
- The device is intended for single use only. Do not resterilize or reuse.
 - Resterilization and/or reuse may result in ineffective catheter coating lubrication, which may result in high friction and the inability to access the target vasculature location; and/or may compromise the structural integrity of the device.
 - Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer/distributor.
 - Use prior to the "Use By" date.
 - Use the BENCHMARK Intracranial Access System in conjunction with fluoroscopic visualization.
 - Do not advance or withdraw the BENCHMARK Intracranial Access System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device. Unrestrained moving or torquing the device against resistance may result in damage to the vessel or device.
 - Maintain a constant infusion of an appropriate flush solution.
 - If flow through the device becomes restricted, do not attempt to clear the lumen by infusion. Remove and replace the device.

Potential adverse events
Possible complications include, but are not limited to, the following: acute occlusion; air embolism; death; distal embolization; emboli; false aneurysm formation; hematoma or hemorrhage at puncture site; infection; intracranial hemorrhage; ischemia; neurological deficits including stroke; vessel spasm, thrombosis, dissection, or perforation.

Penumbra 

www.penumbrainc.com

Penumbra, Inc. USA
One Penumbra Place
Alameda, CA 94502
USA
1.888.272.4606
T 1.510.748.3200
F 1.510.748.3232
order@penumbrainc.com
info@penumbrainc.com

Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician. Prior to use, please refer to the Instructions for Use for complete product indications, contraindications, warnings, precautions, potential adverse events, and detailed instructions for use. Please contact your local Penumbra representative for more information.

Copyright ©2016–2019 Penumbra, Inc. All rights reserved. The Penumbra P logo, Select, Neuron, MAX, and BENCHMARK are registered trademarks or trademarks of Penumbra, Inc. in the USA and other countries. All other trademarks are the property of their respective owners. 11229, Rev. B 03/19 USA