# **COMPASS** Trial

# Aspiration Thrombectomy vs. Stent Retriever Thrombectomy as First-Line Approach

Turk AS, Siddiqui A, Fifi JT, et al. Aspiration thrombectomy versus stent retriever thrombectomy as first-line approach for large vessel occlusion (COMPASS): a multicentre, randomised, open label, blinded outcome, non-inferiority trial. Lancet. 2019;393:998-1008.

Note: COMPASS is an independent, physician-initiated study funded with a research grant by Penumbra, Inc. Penumbra did not play a role in the execution, data collection, data analysis, interpretation, or presentation of results of the COMPASS study.





	ADAPT	SR	OR (95% CI)
Mortality at 90 days	22.0%	22.0%	1.02 (.57, 1.81)
sICH (all ICH with NIHSS ≥ 4 worsening)	6%	6%	1.01 (.37, 2.77)

The opinions and clinical experiences presented herein are for informational purposes only. The results may not be predictive for all patients, Individual results may vary depending on a variety of patient-specific attributes.



\*Per analysis using aggregate supply chain data as primary source and list price as secondary source

# **Trial Background**

# Design

• Prospective, randomized, international, multi-center, blinded assessment concurrent controlled trial

# Population

 Anterior circulation ELVO (ICA to MCA Bifurcation) within 6 hours of onset

# Randomization

1:1 ADAPT vs. Stent Retriever Frontline (SRFL)

# Assessment

- · Blinded core lab adjudication of imaging
- · Blinded mRS certified clinical assessment

# Sites

- 15 centers in North America
  - 6 centers with  $\geq$  67% ADAPT cases
  - 6 centers with  $\geq$  67% SRFL cases
  - 3 centers with mixed technique
  - 270 patients enrolled

# **Baseline Characteristics**

	<b>ADAPT</b> (134)	<b>SRFL</b> (136)	
Age	71.8 ± 13.1	71.1 ± 12.9	
Gender (female)	58%	50%	
Baseline NIHSS (median)	17	17	
Baseline ASPECTS Score (median)	8	8	
Site of Occlusion			
MCA			
M1 Proximal	61% (82/134)	63% (85/136)	
M1 Distal	14% (19/134)	11% (15/136)	
M2 Proximal	8% (11/134)	8% (11/136)	
M3	0 (0/134)	<1% (1/136)	
ICA			
Supraclinoid ICA (ICA Terminus)	13% (18/134)	15% (21/136)	
Petrocavernous	<1% (1/134)	<1% (1/136)	
Other			
Mid-basilar	0 (0/134)	<1% (1/136)	
Tandem Cervical-ICA	2% (3/134)	<1% (1/136)	
Directly admitted to a comprehensive stroke center	56% (75/134)	57% (78/136)	
IV tPA pre-procedure	69% (92/134)	71% (96/136)	
General Anesthesia %	29% (39/134)	30% (41/136)	

Penumbra System - Indication for Use

Penumbra System – Indication for Use Penumbra Reperfusion Catheters and Separators As part of the Penumbra System, the Reperfusion Catheters and Separators are indicated for use in the revascularization of patients with acute ischemic stroke secondary to intracranial large vessel occlusive disease (within the internal carotid, middle cerebral – M1 and M2 segments, basiliar, and vertebral arteries) within B hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment. Penumbra 3D Revascularization Device As part of the Penumbra System, the Penumbra 3D Revascularization Device is indicated for use in the revascularization of patients with acute ischemic stroke secondary to intracranial large vessel occlusive disease (within the internal carotid, middle cerebral – M1 and M2 segments), within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment. Penumbra Aspiration Tubing As part of the Penumbra System, the Penumbra Sterile Aspiration Tubing is indicated to connect the Penumbra Reperfusion

As part of the Penumbra System, the Penumbra Sterile Aspiration Tubing is indicated to connect the Penumbra Reperfusion Catheters to the Penumbra Aspiration Pump. Penumbra Aspiration Pump The Penumbra Aspiration Pump

The Penumbra Aspiration Pump is indicated as a vacuum source for Penumbra Aspiration Systems.

## Contraindications

There are no known contraindications. Warnings

 The Penumbra System should only be used by physicians who have received appropriate training in interventional neuro-endovascular techniques and treatment of acute ischemic stroke.
 Do not advance, retract or use any component of the Penumbra System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device or system as a unit. Unrestrained torquing or forced insertion of the catheter, revascularization device, or separator against resistance may result in damage to the device or

Do not use the Penumbra System with a pump other than the Penumbra Aspiration Pump.

The Penumbra 3D Revascularization Device has not been evaluated in patients with angiographic evidence of pre-existing arterial injury. The use of fluoroscopy may present potential risks from radiation exposure. The probability of adverse events due to radiation The use of fluoroscopy may present potential risks from radiation exposure. The probability of adverse events due to radiation

exposure increases with the total amount of radiation observed, the number of procedures and total procedure time Precautions

- The device is intended for single use only. Do not resterilize or reuse. Resterilization and/or Reuse may result in ineffective catheter coating lubrication, which may result in high friction and the inability to access the target neuro vasculature location.
   Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacture distributive result.
- Use prior to the "Use By" date.
   Use the Penumbra System in conjunction with fluoroscopic visualization.

As in all flournes opported in computation with indexeque precautions to limit patient radiation exposure by using sufficient shielding, reducing fluoroscopy times and modifying radiation technical factors whenever possible.
 Maintain a constant infusion of appropriate flush solution.
 When performing aspiration, ensure that the Penumbra Aspiration Tubing valve is open for only the minimum time needed to

remove thrombus. Excessive aspiration or failure to close the Penumbra Aspiration Tubing valve when aspiration is complete is The Penumbra Separator is not intended for use as a neurovascular guidewire. If repositioning of the Penumbra Reperfusion

Catheter is necessary during the revascularization procedure, such repositioning should be performed over an appropriate neurovascular guidewire using standard microcatheter and guidewire techniques.



Do not use automated high-pressure contrast injection equipment with the Penumbra Reperfusion Catheter because it may damage the device.

- camage the device.
  Administration of anticoagulants and antiplatelets should be suspended until 24 hours posttreatment. Medical management and acute post stroke care should follow the ASA guidelines<sup>1</sup>. Any neurological deterioration should be evaluated by urgent CT scan and other evaluations as indicated according to investigator/hospital best practice.
  As in all surgical interventions, monitoring of intra-procedural blood loss is recommended so that appropriate management may be instituted.
  Confirm vessel diameter, and select an appropriate size Penumbra Reperfusion Catheter. Limit the usage of Reperfusion Catheters to arteries larger than the catheter's outer diameter. Refer to the Reperfusion Catheter labeling for dimensional information.
- information.

Potential Adverse Events Possible complications include, but are not limited to, the following:

Indergic reaction and anaphysication to intrammedia; acute occlusion; air embolism; arteriovenous fistula; death; device mal-function; distal embolization; emboli; false aneurysm formation; hematoma or hemorrhage at access site; inability to completely remove thrombus; infection; intracranal hemorrhage; ischemia; kidney damage from contrast media; neurological deficits including stroke; vessel spasm, thrombosis, dissection, or perforation; radiation exposure that may lead to cataracts, skin reddening or burns from x-ray exposure.

1. Adams, et al., Guidelines for the Early Management of Adults with Ischemic Stroke: A Guideline from the AHA/ASA Stroke Council, Clinical Cardiology Council, Cardiovascular Radiology and Intervention Council, and the Athenosclerotic Peripheral Vascular Disease and Quality of Care Outcomes in Research Interdisciplinary Working Groups: The American Academy of Neurology affirms the value of this guideline as an educational tool for neurologists, Stroke May 2007; 38:1655-1711.

Penumbra Pump MAX – Indication For Use The Penumbra Pump MAX is indicated as a vacuum source for the Penumbra Aspiration Systems. Contraindications

## There are no contraindications.

- Warnings/Precautions
   The canister/tubing is intended for single use only. Do not reuse. Reuse may result in canister cracking or tubing blockages, which may result in the inability to aspirate Do not block bottom or back air vents. Unit may overheat and shut off or fail to restart if run for extended periods of time

- Do not block bottom or back air vents. Unit may overheat and shut off or fail to restart if run for extended periods of time without airflow.
  To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
  Do not position the pump so that it is difficult to operate the power cord disconnection device.
  Remove and service the pump mj flugids or solids have been drawn into the vacuum pump.
  Do not use in the presence of flammable anaesthetic mixture with air or nitrous oxide.
  Do not use in oxygen rich environment.
  To prevent fire or shock hazard, use replacement fuses of equal size and rating.
  To not use pitrole back hazard, use a replacement power cord of equal rating.
  Do not use pitrole back hazard, scies, caustics, or chlorinated solvents to clean or lubricate any parts. It will reduce service life of the pump. Use only water-hase solvents for cleanion. Federal (USA) law restricts this device to sale by or on the order of a physician.
   No modification of this equipment is allowed.

Penumbra, Inc. USA One Penumbra Place Alameda, CA 94502 USA 1.888.272.4606 T 1.510.748.3200 F 1.510.748.3232 order@penumbrainc.com

info@nenumbrainc.com

Penumbra Europe GmbH Am Borsigturm 44 13507 Berlin Germany T +49 30 2005 676-0 F +49 30 2005 676-10 de-order@penumbrainc.com info@penumbrainc.de ıbra Neuro Australia Ptv Ltd

Suite 3, Level 5, 1 Oxford Street Darlinghurst NSW 2010 Australia T +61-1300 817 025 F+61-1300 817 026 order.anz@penumbrainc.com

Product availability varies by country. Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician. Prior to use, please refer to the Instructions for Use for complete product indications, contraindications, warnings, precautions, potential adverse events, and detailed instructions for use. Please contact your local Penumbra representative for more information.

Copyright ©2018–2019 Penumbra, Inc. All rights reserved. The Penumbra P logo, Penumbra System, 3D, 3D Revascularization Device, and MAX are registered trademarks or trademarks of Penumbra, Inc. in the USA and other countries. All other trademarks are the property of their respective owners. 13541, Rev. E 03/19 UNI