


Penumbra 

Artemis™

Neuro Evacuation Device



Stereotactic ICH Underwater Blood Aspiration (SCUBA) technique for minimally invasive endoscopic intracerebral hemorrhage evacuation^a

88.2%
clot evacuation

N = 47 | SD = 20.8

Current Review of ICH Studies: Odds Ratio^b

Comparison	Functional Independence	Survival at Follow-up
MIS vs. Non-MIS	2.2x	1.7x
MIS vs. Craniotomy	2.3x	1.8x
Endoscopic MIS vs. Other Treatments	2.5x	2.7x
Stereotactic Thrombolysis vs. Other Treatment	2.1x	No benefit

a. Kellner CP, Chartrain AG, Nistal DA, et al. The Stereotactic Intracerebral Hemorrhage Underwater Blood Aspiration (SCUBA) technique for minimally invasive endoscopic intracerebral hemorrhage evacuation. J Neurointerv Surg. 2018;10:771-776.
b. Scaggiante J, Zhang X, Mocco J, Kellner CP. Minimally invasive surgery for intracerebral hemorrhage: an updated meta-analysis of randomized controlled trials. Stroke. 2018;49:2612-2620.
The clinical results presented herein are for informational purposes only, and may not be predictive for all patients. Individual results may vary depending on a variety of patient-specific attributes.

Artemis Neuro Evacuation Device



Single Touch Control

Powerful and controlled aspiration

Minimally Invasive Cranial Access

14 mm burr hole with 19 F sheath designed to minimize iatrogenic injury

Neuroendoscope Compatible

Three sizes to work with neuroendoscopes with working channels 1.6–2.9 mm

Atraumatic Evacuation

Recessed bident at distal tip maintains cannula patency and aspiration


Atraumatic Bident

Controlled Aspiration


Maximum vacuum power
-29 inHg^a



Sheath Insertion





Active Evacuation





Cavity Exploration

a. Tests performed and data on file at Penumbra, Inc. Bench test results may not be indicative of clinical performance.


Case Examples





Christopher Kellner, MD
Mount Sinai Health, NY




Dimitri Sigounas, MD
George Washington University Medical Center, Washington, D.C.





Paul Saphier, MD
Overlook Medical Center, NJ



Ziad Hage, MD
Novant Health, NC



Christopher Nickle, MD
Methodist University Hospital, TN



Pinakin R. Jethwa, MD
Atlantic NeuroSurgical Specialists, NJ

Images used with permission. Consents on file at Penumbra, Inc. Images for illustrative purposes only. Individual results may vary depending on a variety of patient-specific attributes.



Artemis™

Neuro Evacuation Device

Ordering Information

Catalog Number	Description	Cannula Outer Diameter	Cannula Length	Aspiration Tubing Length
AP28	Artemis Device 2.8 mm	2.8 mm	27 cm	9.5 ft
AP21	Artemis Device 2.1 mm	2.1 mm	26 cm	9.5 ft
AP15	Artemis Device 1.5 mm	1.5 mm	27 cm	9.5 ft
APCAN2	MAX Canister	—	—	—
PMX110	Pump MAX 110V	—	—	—

ARTEMIS Neuro Evacuation Device — Indication for Use

The ARTEMIS Neuro Evacuation Device is used for the controlled aspiration of tissue and/or fluid during surgery of the Ventricular System or Cerebrum in conjunction with a Penumbra Aspiration Pump.

Penumbra Aspiration Pump:

The Penumbra Aspiration Pump is indicated as a vacuum source for the Penumbra Aspiration Systems.

Contraindications

- The ARTEMIS Neuro Evacuation Device is not recommended during surgery of the brainstem, cerebellum, epidural or subdural spaces. • Do not use fibrinolytic therapy during the procedure. • Do not use the ARTEMIS Neuro Evacuation Device with a non-Penumbra recommended Aspiration Pump. The safety and effectiveness of its use with a non-Penumbra recommended Aspiration Pump has not been established and can lead to patient injury or death.

Warnings

The ARTEMIS Neuro Evacuation Device should only be used by physicians who have received appropriate training to perform image-guided neurosurgical procedures.

Precautions

- The ARTEMIS Neuro Evacuation Device is intended for single use only. Do not resterilize or reuse. Resterilization or reuse could lead to infection or ineffective

removal of tissue and/or fluid. • Do not use kinked or damaged devices. Do not use open or damaged packages. Return damaged devices and packaging to the manufacturer/distributor. • Use prior to the "Use By" date. • Use the ARTEMIS Neuro Evacuation Device in conjunction with intraprocedural image-guidance. • Do not use in an oxygen rich environment. • Do not advance or use the ARTEMIS Neuro Evacuation Device against resistance without careful visual assessment of the cause. If the cause cannot be determined, withdraw the device. Unrestrained torqueing or forced insertion of the device against resistance may result in damage to the device, which may lead to tissue damage and/or device breakage.

Potential Adverse Events

Possible complications include, but are not limited to, the following: hematoma expansion, fever, headaches, vomiting, hyperglycemia, edema, re-bleeding, death, bleeding, increased blood pressure, infections, seizures, intraventricular hemorrhage, hydrocephalus, thromboembolic events, decreased consciousness, craniotomy, unintended removal of tissue leading to neurological and/or sensory deficit.

Penumbra Pump MAX — Indication for Use

The Penumbra Pump MAX is indicated as a vacuum source for the Penumbra Aspiration Systems.

Contraindications

There are no contraindications.

Warnings/Precautions

- The canister/tubing is intended for single use only. Do not reuse. Reuse may result in canister cracking or tubing blockages, which may result in the inability to aspirate. • Do not block bottom or back air vents. Unit may overheat and shut off or fail to restart if run for extended periods of time without airflow. • To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth. • Do not position the pump so that it is difficult to operate the power cord disconnection device. • Remove and service the pump if liquids or solids have been drawn into the vacuum pump. • Do not use in the presence of flammable anaesthetic mixture with air or nitrous oxide. • Do not use in oxygen rich environment. • To prevent fire or shock hazard, use replacement fuses of equal size and rating. • To prevent fire or shock hazard, use a replacement power cord of equal rating. • Do not re-infuse blood or fluid from the canister back into the patient. • Do not use petroleum base compounds, acids, caustics, or chlorinated solvents to clean or lubricate any parts. It will reduce service life of the pump. Use only water-base solvents for cleaning. • Federal (USA) law restricts this device to sale by or on the order of a physician. • No modification of this equipment is allowed.

Penumbra

www.penumbrainc.com

Penumbra, Inc. USA
One Penumbra Place
Alameda, CA 94502
USA
1.888.272.4606
T 1.510.748.3200
F 1.510.748.3232
order@penumbrainc.com
info@penumbrainc.com

Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician. Prior to use, please refer to the Instructions for Use for complete product indications, contraindications, warnings, precautions, potential adverse events, and detailed instructions for use. Rendering for illustrative purposes only. Individual results may vary depending on a variety of patient-specific attributes. Please contact your local Penumbra representative for more information.

Copyright ©2017–2019 Penumbra, Inc. All rights reserved. The Penumbra P logos, Artemis, and MAX are registered trademarks or trademarks of Penumbra, Inc. in the USA and other countries. 12562, Rev. D 09/19 USA