

SIM-V

SIM

BER

H1



Seamless transition zone with no shelf for atraumatic advancement




### Available Lengths

<b>5F</b>	–	130 cm	120 cm 130 cm	120 cm 130 cm
<b>6F</b>	125 cm	125 cm	105 cm 125 cm	105 cm 125 cm

## All Penumbra Select Catheters


- Stainless braided shaft with radiopaque polymer steam-shaped tip
- Tapered profile of each Select catheter optimized for seamless transition and torque response
- Penumbra Select catheters do not have hydrophilic coating

### 5F Select | 120/130 cm



- Designed to deliver soft-tipped BENCHMARK® 071 and Neuron™ 070 intracranial access catheters
- Not designed to deliver 6F long sheath

### 6F Select | 105/125 cm



- Designed to deliver Neuron MAX long sheath
- Maximum OD spec = .087"
- Not designed for use with smaller ID catheters