

Penumbra

DDC

Hyper Distal Delivery Catheter

5MAX™ DDC

Compatibility

5MAX DDC 115 is compatible with ALL neurovascular microcatheter, balloon and stent/flow diverter delivery systems.

Pericallosal AVM



Neuron™ MAX 088

Penumbra 

www.penumbrainc.com

Penumbra DDC Hyper Distal Delivery Catheter

Penumbra **Hyper Distal Delivery Catheters** are designed for superior tracking and angiographic runs to simplify therapy delivery to the most distal, tortuous anatomy. The unique Penumbra tapered inner shaft maximizes the effective flow lumen to deliver exceptional roadmapping. The short lengths enable a wide variety of device compatibility.



DDC Ordering Information 5MAX DDC 115 is compatible with ALL neurovascular microcatheter, balloon and stent/FD delivery systems.

Catalog Number	Description	Proximal OD	Distal OD	Effective Flow Lumen	Working Length	Total Length
5MAXDDC115	5MAX™ Distal Delivery Catheter - 115 cm	6.0F (.080")	5.0F	.062"	115 cm	119 cm
5MAXDDC125	5MAX™ Distal Delivery Catheter - 125 cm	6.0F (.080")	5.0F	.062"	125 cm	129 cm
4MAXDDC130	4MAX™ Distal Delivery Catheter - 130 cm	6.0F (.080")	4.3F	.053"	130 cm	134 cm

Neuron™ MAX 6F 088 + 6F Select™ Catheter Recommended

Catalog Number	Description	Tip Shape	Proximal OD ^a	Distal OD ^a	Working Length	Inner Diameter
Neuron™ MAX 6F 088 Lumen Long Sheath (Crosscut Valve, RHV and Dilator included)						
PNML6F088804	6F 088 Neuron™ MAX Long Sheath, 80/4 Straight		8.0F	8.0F	80 cm	.088"
PNML6F088804M	6F 088 Neuron™ MAX Long Sheath, 80/4 MP		8.0F	8.0F	80 cm	.088"
PNML6F088904	6F 088 Neuron™ MAX Long Sheath, 90/4 Straight		8.0F	8.0F	90 cm	.088"
PNML6F088904M	6F 088 Neuron™ MAX Long Sheath, 90/4 MP		8.0F	8.0F	90 cm	.088"
6F Select™ Catheter						
PNS6F105BER	6F Select™ Catheter, 105 BER		5.6F / 6.0F	5.0F	105 cm	.040"
PNS6F105H1	6F Select Catheter, 105 H1		5.6F / 6.0F	5.0F	105 cm	.040"
PNS6F125BER	6F Select Catheter, 125 BER		5.6F / 6.0F	5.0F	125 cm	.040"
PNS6F125SIM	6F Select Catheter, 125 SIM		5.6F / 6.0F	5.0F	125 cm	.040"
PNS6F125H1	6F Select Catheter, 125 H1		5.6F / 6.0F	5.0F	125 cm	.040"
PNS6F125SIMV	6F Select Catheter, 125 SIM-V		5.6F / 6.0F	5.0F	125 cm	.040"

a. A mid outer diameter is only listed if applicable to device.



www.penumbrainc.com

Penumbra, Inc. USA
 One Penumbra Place
 Alameda, CA 94502
 USA
 1.888.272.4606
 T 1.510.748.3200
 F 1.510.748.3232
 order@penumbrainc.com
 info@penumbrainc.com

Penumbra Europe GmbH
 Am Borsigturm 44
 13507 Berlin
 Germany
 T +49 30 2005 676-0
 F +49 30 2005 676-10
 order@penumbrainc.de
 info@penumbrainc.de

Penumbra Neuro Australia Pty Ltd
 Suite 3, Level 5, 1 Oxford Street
 Darlinghurst NSW 2010
 Australia
 T +61-1300 817 025
 F +61-1300 817 026
 order.anz@penumbrainc.com