

Penumbra System®

ELEVATING PERFORMANCE WITH

Penumbra

JET 7™

Penumbra

+ ENGINE™



Penumbra 

Penumbra
JET 7



.072" LUMEN

20 TRANSITIONS
for trackability and navigation


ARTICULATING MARKER BAND
designed to improve tip softness

SUPERIOR FLEXIBILITY
enabled by progressive distal coil wind

ENHANCED PUSHABILITY
featuring Quad-Wire technology

FULL LENGTH PTFE LINER
designed for durability with adjunctive devices

2.16 mm




132 cm length

HIGHEST TRF

WITH PENUMBRA JET 7
POWERED BY PENUMBRA ENGINE

$$\text{TRF} = \frac{\text{Catheter Tip Area}}{\text{Vacuum Level}} \times$$

Thrombus removal force (TRF) equates to catheter tip force and drives clot capture.



a. Tests performed and data on file at Penumbra, Inc. Bench test results may not be indicative of clinical performance.

